


## How to set the parameters for single point confocals

- 1. Scientific question:** Define the scientific question for this image. Usually 1 question per image. What do you want to measure? Which resolution, working distance, field of view and immersion medium are needed?
- 2. Fluorophores:** Choose your fluorophore combination: match fluorophore spectra with the microscope excitation and emission specifications.
- 3. Choose an imaging strategy:**
  - Imaging sequentially? Simultaneously? Combining channels? Spectrally?
  - Avoid bleedthrough while optimizing the acquisition time.
- 4. Objective:** Choose the objective that delivers the required resolution and fulfils the other criteria ( $0.61 * \lambda / NA$ ).
- 5. Display:**
  - Set the images to be shown in black and white to avoid colour bias.
  - Use the LUT to push the brightness and contrast.
  - Check the pixel intensity in and out of the region of interest to estimate the signal to background ratio.
- 6. Zoom to adjust the scan area:** The area of interest should fill the image, no unnecessary information.
- 7. Pixel density:**
  - What is the pixel size that fulfills the Nyquist criterion? Click on Optimal/Nyquist to find out ( $resolution/2$ ).
  - Do you need to use such small pixels to reliably answer your scientific question? Yes, if small distances must be measured.
  - Can you afford to use such small pixels? Only if bleaching, light toxicity... are not an issue.
  - Beware of aliasing/Moiré artifacts for small details: Try different degrees of undersampling and compare reliably the images answer your scientific question.
- 8. Pinhole size:**
  - Should the pinhole to be optimized for the shorter or longer wavelength?
    - Best contrast for short wavelength -> dimmer for long wavelength.
    - Best contrast for long wavelength -> lower contrast for short wavelength.
  - Open the pinhole  $>1$  AU for very dim samples that do not require high contrast in xyz.
  - Close the pinhole  $<1$  AU if very bright samples to get higher xy and especially z contrast.



## 9. Laser power:

- Stay low if bleaching or light toxicity are a problem (live samples, z stacks, dim samples).
- Increase to reduce shot noise contribution and improve the S/N ratio.
- Do not saturate the fluorophores (it decreases the S/N ratio). Check by doubling laser power. The fluorescence intensity should double.
- Do not send the fluorophores to the dark state (fluorophore depletion which decreases the S/N ratio).

## 10. Detector gain and offset

- High gain creates noise. Set a low to medium gain and compensate by increasing the laser power.
- Set the offset to have no pixel with a 0 intensity value in the background (using the saturation/underexposure colour mask).

## 11. Noise:

- Acquire a timelapse of 3 images without interval/delay. What changes from one image to the next is random noise.
- Does the random noise prevent you from answering your scientific question? If yes, improve the S/N ratio by increasing the illumination intensity, averaging or using a denoising post-acquisition tool.

## 12. Pixel dwell time

- If all other noise-reducing strategies are insufficient, lower the scan speed (= increase the pixel dwell time). More photons go to the detector so less amplification (gain) is needed.
- Disadvantages: more bleaching, longer imaging, higher light toxicity, higher risk of fluorophore depletion (dark state).