

# BD LSR II User Guide

## Starting up

1. **Switch on the Fluidics Cart** (located under the cytometer)
2. Turn on the power to the cytometer
3. Check the fluid levels, that the FACSFlow tank is full, and the Waste tank is empty.  
*\*In case of exchanging any of the tanks, press Restart in the Fluidics Cart right after.*
4. Switch on the computer, log in to Windows under the BDOperator user, launch and log in into the Diva Software. Wait until change to connected status.
5. Pre-run:
  - a. Prime twice without tube on Sample Injection Port (SIP) (pressing PRIME in the control panel of the cytometer).
  - b. Pre-run Cleaning: run FACSRinse (5min) followed by diH<sub>2</sub>O (5min) while you set up for your experiment. Run at a HIGH flow rate.
  - c. In case of being the first user of the day, run CS&T performance check (**Cytometer > CST**). Use an already prepared tube with CST-beads (good for up to 5 days if stored cold) or prepare a new one according to the instructions on the box.
6. The instrument is ready to run your samples.

## Shutting down

1. Copy your files to an external USB (data will be deleted by the BFC every two weeks).
2. **After-run cleaning:** all washes should be done in **HIGH flow rate:**
  - a. Run FACSClean with the lever open until tube is empty
  - b. Run 10min FACSClean with the lever closed
  - c. Run FACSRinse with the lever open until tube is empty
  - d. Run 10min FACSRinse with the lever closed
  - e. Run diH<sub>2</sub>O with the lever open until tube is empty
  - f. Run 10min diH<sub>2</sub>O with the lever closed
  - g. Leave a tube with diH<sub>2</sub>O on the sample probe
  - h. Make sure the instrument is set to Standby mode
3. Exit Diva software
4. Turn off cytometer power
5. Turn off the computer
6. Turn off the Fluidics Cart

## Frequent problems LSR II

Feel free to ask BFC staff if you can not solve your issue using this guide

### **Problem: Alarm sounds on Fluidics Cart**

#### **Causes and solutions:**

- **Full waste container OR Empty sheath container:**
  - Switch off the Fluidics Cart
  - Replace the full waste container by a new empty box OR get a new full BD FACS Sheath Container
  - Restart Fluidics Cart
  - Prime the fluidics after changing the tanks:
    - Install an empty tube on the SIP
    - Press twice the PRIME fluid control button
    - Install a FACS tube with 1ml of deionised water on the SIP and run it on high flow rate for about 2min.
- **Waste OR sheath detectors are dirty**
  - Remove the lid of the Waste OR Sheath container and clean the detectors with a wipe and ethanol.
  - Screw up the lid again
  - Restart the Fluidics Cart

### **Problem: Very low or no acquisition rate**

#### **Causes and solutions:**

- **Cracked FACS tube**
  - Transfer the content to a new FACS tube and rerun the sample.
  - If there is no acquisition, check for blockage or air lock.
- **System blockage**
  - Blockage needs to be shifted:
    - Turn the instrument to Standby and try to unblock the machine by pressing the Prime button.
    - If that did not work, run some hot water.

### **Problem: No events shown in the dotplots**

#### **Causes and solutions:**

- **Diva Software is running slow**
  - Press Ctrl+Shift+Delete to open the Task Manager. Close the Task Manager and go back to Diva.
  - In case events are not shown, restart the computer and try again.