

Request FENO

Study nr: _____

(Filled in by FENO)

Name	Cost site (ZZ-nr. & project nr.)
Institution /Address	Date

Species/ nr. of animals
Ethical permit nr.
Animal ID-nr. & tissues
Examination

Study nr: _____

General instructions for histological sampling (IF using 4% formaldehyde):

- Contact us beforehand.
- General fixation of tissues in 4% neutral buffered formaldehyde or freshly prepared paraformaldehyde (PFA) for a preferred maximum time of 24 h. To ensure optimal fixation the volume of the fixative should be at least ten times the volume of the tissue, and tissue parts should not be thicker than 5 mm. Use flat bottom vials/jars.
- After 24 hours transfer tissues to 70 % ethanol unless working with brain or adipose.

Filled in by the FENO lab:

Slides/blocks delivered to pathologist/customer:

Date

Responsible pathologist/customer

Responsible technician

	Number		Number
VIP		IHC	
Blocks		Special stain	
H&E		Tissue remains	
Plus glass		Archive	

Comments: _____