

Well #	Salt	Well #	Buffer ◇	Well #	Precipitant
1. (A1)	0.1 M Sodium chloride	1. (A1)	0.1 M Sodium acetate trihydrate pH 4.6	1. (A1)	12% v/v (+/-)-2-Methyl-2,4-pentanediol
2. (A2)	0.1 M Zinc acetate dihydrate	2. (A2)	0.1 M Sodium acetate trihydrate pH 4.6	2. (A2)	12% w/v Polyethylene glycol 4,000
3. (A3)	0.2 M Ammonium sulfate	3. (A3)	0.1 M Sodium acetate trihydrate pH 4.6	3. (A3)	10% w/v Polyethylene glycol 4,000
4. (A4)	0.1 M Sodium chloride	4. (A4)	0.1 M Sodium acetate trihydrate pH 4.6	4. (A4)	12% v/v 2-Propanol
5. (A5)	None	5. (A5)	0.1 M Sodium acetate trihydrate pH 4.6	5. (A5)	12% w/v Polyethylene glycol 4,000
6. (A6)	None	6. (A6)	0.1 M Sodium acetate trihydrate pH 4.6	6. (A6)	1.0 M Ammonium sulfate
7. (A7)	None	7. (A7)	0.1 M Sodium acetate trihydrate pH 4.6	7. (A7)	1.0 M Magnesium sulfate heptahydrate
8. (A8)	0.1 M Magnesium chloride hexahydrate	8. (A8)	0.1 M Sodium acetate trihydrate pH 4.6	8. (A8)	18% v/v Polyethylene glycol 400
9. (A9)	0.1 M Lithium sulfate monohydrate	9. (A9)	0.1 M Sodium acetate trihydrate pH 4.6	9. (A9)	1.0 M Ammonium phosphate monobasic
10. (A10)	0.1 M Sodium chloride	10. (A10)	0.1 M Sodium acetate trihydrate pH 4.6	10. (A10)	12% w/v Polyethylene glycol 6,000
11. (A11)	0.1 M Magnesium chloride hexahydrate	11. (A11)	0.1 M Sodium acetate trihydrate pH 4.6	11. (A11)	12% w/v Polyethylene glycol 6,000
12. (A12)	0.1 M Sodium chloride	12. (A12)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	12. (A12)	18% v/v Polyethylene glycol 400
13. (B1)	0.1 M Lithium sulfate monohydrate	13. (B1)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	13. (B1)	12% w/v Polyethylene glycol 4,000
14. (B2)	0.1 M Sodium citrate tribasic dihydrate	14. (B2)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	14. (B2)	10% v/v 2-Propanol
15. (B3)	0.1 M Sodium chloride	15. (B3)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	15. (B3)	12% v/v (+/-)-2-Methyl-2,4-pentanediol
16. (B4)	None	16. (B4)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	16. (B4)	1.0 M Magnesium sulfate heptahydrate
17. (B5)	0.1 M Sodium chloride	17. (B5)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	17. (B5)	12% w/v Polyethylene glycol 4,000
18. (B6)	0.1 M Lithium sulfate monohydrate	18. (B6)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	18. (B6)	12% w/v Polyethylene glycol 6,000
19. (B7)	0.1 M Magnesium chloride hexahydrate	19. (B7)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	19. (B7)	4% v/v (+/-)-2-Methyl-2,4-pentanediol
20. (B8)	None	20. (B8)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	20. (B8)	0.1 M Sodium chloride
21. (B9)	0.1 M Lithium sulfate monohydrate	21. (B9)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	21. (B9)	4% v/v Polyethylene glycol 400
22. (B10)	None	22. (B10)	0.1 M ADA pH 6.5	22. (B10)	1.0 M Ammonium sulfate
23. (B11)	0.1 M Lithium sulfate monohydrate	23. (B11)	0.1 M ADA pH 6.5	23. (B11)	12% w/v Polyethylene glycol 4,000, 2% v/v 2-Propanol
24. (B12)	None	24. (B12)	0.1 M ADA pH 6.5	24. (B12)	1.0 M Ammonium phosphate dibasic
25. (C1)	0.1 M Magnesium chloride hexahydrate	25. (C1)	0.1 M ADA pH 6.5	25. (C1)	12% w/v Polyethylene glycol 6,000
26. (C2)	None	26. (C2)	0.1 M ADA pH 6.5	26. (C2)	12% v/v (+/-)-2-Methyl-2,4-pentanediol
27. (C3)	0.1 M Lithium sulfate monohydrate	27. (C3)	0.1 M ADA pH 6.5	27. (C3)	1.0 M Magnesium sulfate hydrate
28. (C4)	0.3 M Lithium sulfate monohydrate	28. (C4)	0.1 M ADA pH 6.5	28. (C4)	4% v/v Polyethylene glycol 400
29. (C5)	0.1 M Ammonium sulfate	29. (C5)	0.1 M HEPES sodium pH 7.5	29. (C5)	0.5 M Sodium phosphate dibasic dihydrate, 0.5 M Potassium phosphate dibasic
30. (C6)	0.1 M Sodium chloride	30. (C6)	0.1 M HEPES sodium pH 7.5	30. (C6)	10% w/v Polyethylene glycol 4,000
31. (C7)	0.1 M Magnesium chloride hexahydrate	31. (C7)	0.1 M HEPES sodium pH 7.5	31. (C7)	18% v/v Polyethylene glycol 400
32. (C8)	None	32. (C8)	0.1 M HEPES sodium pH 7.5	32. (C8)	1.0 M Potassium sodium tartrate tetrahydrate
33. (C9)	0.1 M Ammonium sulfate	33. (C9)	0.1 M HEPES sodium pH 7.5	33. (C9)	18% v/v Polyethylene glycol 400
34. (C10)	0.1 M Ammonium sulfate	34. (C10)	0.1 M HEPES sodium pH 7.5	34. (C10)	10% w/v Polyethylene glycol 4,000
35. (C11)	0.1 M Sodium citrate tribasic dihydrate	35. (C11)	0.1 M HEPES sodium pH 7.5	35. (C11)	12% v/v (+/-)-2-Methyl-2,4-pentanediol
36. (C12)	None	36. (C12)	0.1 M HEPES sodium pH 7.5	36. (C12)	1.0 M Sodium citrate tribasic dihydrate
37. (D1)	0.6 M Magnesium sulfate hydrate	37. (D1)	0.1 M HEPES sodium pH 7.5	37. (D1)	4% v/v Polyethylene glycol 400
38. (D2)	0.6 M Magnesium sulfate hydrate	38. (D2)	0.1 M HEPES sodium pH 7.5	38. (D2)	4% v/v (+/-)-2-Methyl-2,4-pentanediol
39. (D3)	0.1 M Lithium sulfate monohydrate	39. (D3)	0.1 M HEPES sodium pH 7.5	39. (D3)	0.1 M Potassium sodium tartrate tetrahydrate
40. (D4)	0.1 M Lithium sulfate monohydrate	40. (D4)	0.1 M TRIS hydrochloride pH 8.5	40. (D4)	12% v/v (+/-)-2-Methyl-2,4-pentanediol
41. (D5)	0.1 M Ammonium phosphate dibasic	41. (D5)	0.1 M TRIS hydrochloride pH 8.5	41. (D5)	0.5 M Sodium phosphate dibasic dihydrate, 0.5 M Potassium phosphate dibasic
42. (D6)	None	42. (D6)	0.1 M TRIS hydrochloride pH 8.5	42. (D6)	0.1 M Sodium acetate trihydrate
43. (D7)	None	43. (D7)	0.1 M TRIS hydrochloride pH 8.5	43. (D7)	0.1 M Sodium chloride
44. (D8)	0.1 M Ammonium phosphate dibasic	44. (D8)	0.1 M TRIS hydrochloride pH 8.5	44. (D8)	12% w/v Polyethylene glycol 6,000
45. (D9)	0.1 M Potassium sodium tartrate tetrahydrate	45. (D9)	0.1 M TRIS hydrochloride pH 8.5	45. (D9)	0.4 M Magnesium sulfate hydrate
46. (D10)	None	46. (D10)	0.1 M TRIS hydrochloride pH 8.5	46. (D10)	0.2 M Lithium sulfate monohydrate
47. (D11)	None	47. (D11)	0.1 M TRIS hydrochloride pH 8.5	47. (D11)	0.5 M Ammonium sulfate
48. (D12)	0.1 M Sodium citrate tribasic dihydrate	48. (D12)	0.1 M TRIS hydrochloride pH 8.5	48. (D12)	5% v/v Polyethylene glycol 400

◇ Buffer pH is that of a 1.0 M stock prior to dilution with other reagent components: pH with HCl or NaOH.

*MembFac HT™ (Deep Well Block) contains ninety-six unique reagents beginning at position A1.  
To determine the formulation of each reagent, simply read across the page.*

Well #	Salt	Well #	Buffer ◇	Well #	Precipitant
49. (E1)	0.02 M Calcium chloride dihydrate	49. (E1)	0.1 M Sodium acetate trihydrate pH 4.6	49. (E1)	15% v/v (+/-)-2-Methyl-2,4-pentanediol
50. (E2)	None	50. (E2)	None	50. (E2)	0.2 M Potassium sodium tartrate tetrahydrate
51. (E3)	None	51. (E3)	None	51. (E3)	0.2 M Ammonium phosphate monobasic
52. (E4)	None	52. (E4)	0.1 M TRIS hydrochloride pH 8.5	52. (E4)	1.0 M Ammonium sulfate
53. (E5)	0.2 M Sodium citrate tribasic dihydrate	53. (E5)	0.1 M HEPES sodium pH 7.5	53. (E5)	15% v/v (+/-)-2-Methyl-2,4-pentanediol
54. (E6)	0.2 M Magnesium chloride hexahydrate	54. (E6)	0.1 M TRIS hydrochloride pH 8.5	54. (E6)	15% w/v Polyethylene glycol 4,000
55. (E7)	None	55. (E7)	0.1 M Sodium cacodylate trihydrate pH 6.5	55. (E7)	0.7 M Sodium acetate trihydrate
56. (E8)	0.2 M Sodium citrate tribasic dihydrate	56. (E8)	0.1 M Sodium cacodylate trihydrate pH 6.5	56. (E8)	15% v/v 2-Propanol
57. (E9)	0.2 M Ammonium acetate	57. (E9)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	57. (E9)	15% w/v Polyethylene glycol 4,000
58. (E10)	0.2 M Ammonium acetate	58. (E10)	0.1 M Sodium acetate trihydrate pH 4.6	58. (E10)	15% w/v Polyethylene glycol 4,000
59. (E11)	None	59. (E11)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	59. (E11)	0.5 M Ammonium phosphate monobasic
60. (E12)	0.2 M Magnesium chloride hexahydrate	60. (E12)	0.1 M HEPES sodium pH 7.5	60. (E12)	15% v/v 2-Propanol
61. (F1)	0.2 M Sodium citrate tribasic dihydrate	61. (F1)	0.1 M TRIS hydrochloride pH 8.5	61. (F1)	15% v/v Polyethylene glycol 400
62. (F2)	0.2 M Calcium chloride dihydrate	62. (F2)	0.1 M HEPES sodium pH 7.5	62. (F2)	14% v/v Polyethylene glycol 400
63. (F3)	0.2 M Ammonium sulfate	63. (F3)	0.1 M Sodium cacodylate trihydrate pH 6.5	63. (F3)	15% w/v Polyethylene glycol 8,000
64. (F4)	None	64. (F4)	0.1 M HEPES sodium pH 7.5	64. (F4)	0.75 M Lithium sulfate monohydrate
65. (F5)	0.2 M Lithium sulfate monohydrate	65. (F5)	0.1 M TRIS hydrochloride pH 8.5	65. (F5)	15% w/v Polyethylene glycol 4,000
66. (F6)	0.2 M Magnesium acetate tetrahydrate	66. (F6)	0.1 M Sodium cacodylate trihydrate pH 6.5	66. (F6)	10% w/v Polyethylene glycol 8,000
67. (F7)	0.2 M Ammonium acetate	67. (F7)	0.1 M TRIS hydrochloride pH 8.5	67. (F7)	15% v/v 2-Propanol
68. (F8)	0.2 M Ammonium sulfate	68. (F8)	0.1 M Sodium acetate trihydrate pH 4.6	68. (F8)	12.5% w/v Polyethylene glycol 4,000
69. (F9)	0.2 M Magnesium acetate tetrahydrate	69. (F9)	0.1 M Sodium cacodylate trihydrate pH 6.5	69. (F9)	15% v/v (+/-)-2-Methyl-2,4-pentanediol
70. (F10)	0.2 M Sodium acetate trihydrate	70. (F10)	0.1 M TRIS hydrochloride pH 8.5	70. (F10)	15% w/v Polyethylene glycol 4,000
71. (F11)	0.2 M Magnesium chloride hexahydrate	71. (F11)	0.1 M HEPES sodium pH 7.5	71. (F11)	15% v/v Polyethylene glycol 400
72. (F12)	0.2 M Calcium chloride dihydrate	72. (F12)	0.1 M Sodium acetate trihydrate pH 4.6	72. (F12)	10% v/v 2-Propanol
73. (G1)	None	73. (G1)	0.1 M Imidazole pH 6.5	73. (G1)	0.5 M Sodium acetate trihydrate
74. (G2)	0.2 M Ammonium acetate	74. (G2)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	74. (G2)	15% v/v (+/-)-2-Methyl-2,4-pentanediol
75. (G3)	0.2 M Sodium citrate tribasic dihydrate	75. (G3)	0.1 M HEPES sodium pH 7.5	75. (G3)	10% v/v 2-Propanol
76. (G4)	0.2 M Sodium acetate trihydrate	76. (G4)	0.1 M Sodium cacodylate trihydrate pH 6.5	76. (G4)	15% w/v Polyethylene glycol 8,000
77. (G5)	None	77. (G5)	0.1 M HEPES sodium pH 7.5	77. (G5)	0.4 M Potassium sodium tartrate tetrahydrate
78. (G6)	0.2 M Ammonium sulfate	78. (G6)	None	78. (G6)	15% w/v Polyethylene glycol 8,000
79. (G7)	0.2 M Ammonium sulfate	79. (G7)	None	79. (G7)	15% w/v Polyethylene glycol 4,000
80. (G8)	None	80. (G8)	None	80. (G8)	1.0 M Ammonium sulfate
81. (G9)	None	81. (G9)	None	81. (G9)	2.0 M Sodium formate
82. (G10)	None	82. (G10)	0.1 M Sodium acetate trihydrate pH 4.6	82. (G10)	1.0 M Sodium formate
83. (G11)	None	83. (G11)	0.1 M HEPES sodium pH 7.5	83. (G11)	0.4 M Sodium phosphate monobasic monohydrate, 0.4 M Potassium phosphate monobasic
84. (G12)	None	84. (G12)	0.1 M TRIS hydrochloride pH 8.5	84. (G12)	4% w/v Polyethylene glycol 8,000
85. (H1)	None	85. (H1)	0.1 M Sodium acetate trihydrate pH 4.6	85. (H1)	4% w/v Polyethylene glycol 4,000
86. (H2)	None	86. (H2)	0.1 M HEPES sodium pH 7.5	86. (H2)	0.7 M Sodium citrate tribasic dihydrate
87. (H3)	None	87. (H3)	0.1 M HEPES sodium pH 7.5	87. (H3)	2% v/v Polyethylene glycol 400, 1.0 M Ammonium sulfate
88. (H4)	None	88. (H4)	0.1 M Sodium citrate tribasic dihydrate pH 5.6	88. (H4)	10% v/v 2-Propanol, 10% w/v Polyethylene glycol 4,000
89. (H5)	None	89. (H5)	0.1 M HEPES sodium pH 7.5	89. (H5)	5% v/v 2-Propanol, 10% w/v Polyethylene glycol 4,000
90. (H6)	0.05 M Potassium phosphate monobasic	90. (H6)	None	90. (H6)	10% w/v Polyethylene glycol 8,000
91. (H7)	None	91. (H7)	None	91. (H7)	15% w/v Polyethylene glycol 1,500
92. (H8)	None	92. (H8)	None	92. (H8)	0.1 M Magnesium formate dihydrate
93. (H9)	0.2 M Zinc acetate dihydrate	93. (H9)	0.1 M Sodium cacodylate trihydrate pH 6.5	93. (H9)	9% w/v Polyethylene glycol 8,000
94. (H10)	0.2 M Calcium acetate hydrate	94. (H10)	0.1 M Sodium cacodylate trihydrate pH 6.5	94. (H10)	9% w/v Polyethylene glycol 8,000
95. (H11)	None	95. (H11)	0.1 M Sodium acetate trihydrate pH 4.6	95. (H11)	1.0 M Ammonium sulfate
96. (H12)	None	96. (H12)	0.1 M TRIS hydrochloride pH 8.5	96. (H12)	1.0 M Ammonium phosphate monobasic

◇ Buffer pH is that of a 1.0 M stock prior to dilution with other reagent components: pH with HCl or NaOH.

*MembFac HT™ (Deep Well Block) contains ninety-six unique reagents beginning at position A1.  
To determine the formulation of each reagent, simply read across the page.*

**HAMPTON  
RESEARCH**

*Solutions for Crystal Growth*

34 Journey

Aliso Viejo, CA 92656-3317 U.S.A.

Tel: (949) 425-1321 • Fax: (949) 425-1611

E-mail: tech@hrmail.com

Website: www.hamptonresearch.com

© 1991-2016 Hampton Research Corp. all rights reserved  
Printed in the United States of America. This guide or  
parts thereof may not be reproduced in any form without  
the written permission of the publishers.