

BD FACSCanto™ II

Laser	Filter	LP	PMT	Examples:
Violet : 405	B	450/50	DAPI	Pacific blue, BV421, V450, Violet Live/Dead, SBV440, VioBrightV423, VioBlue, eF450
	A	510/50	505LP AmCyan	BV510, V500, Aqua Live/Dead, SBV475, SBV515, BioGreen, ZombieAqua
Blue: 488	E	488/10	SSC	
	D	530/30	502LP FITC	Alexa488, GFP, Green Live/Dead, BB515, VioBright515, VioBrightFITC, CFSE
	C	585/42	556LP PE	tdTomato, dsRed, SBY575
	B	670LP	655LP PE-Cy5	7AAD, PE-A647, SBY665
	A	780/60	735LP PE-Cy7	PE-A750, SBY800, PE-Vio 770
Red: 640	B	660/42	Alexa 647	Alexa 647, Far Red Live/Dead, SBR670, eF660, Vio Bright R667
	A	780/60	735LP APC-Cy7	APC-H7, APC-eF780, SBR775, APC-Fire750, APC-Vio 770, L/D NIR, Zombie NIR