

BD FACSVeTM

Laser		Filter	LP	PMT	Examples:
Violet : 405	B	448/45		DAPI	Pacific blue, BV421, V450, Violet Live/Dead, SBV440, VioBrightV423, VioBlue, eF450
	A	528/45		AmCyan	BV510, V500, Aqua Live/Dead, SBV475, SBV515, BioGreen, ZombieAqua
Blue: 488	E	488/15		SSC	
	D	527/32		FITC	Alexa488, GFP, Green Live/Dead, BB515, VioBright515, VioBrightFITC, CFSE
	C	586/42		PE	tdTomato, dsRed, SBY575
	B	700/54		PE-Cy5	PerCP-Cy5.5, PerCP, PerCP-eF710, SBB675, SBB700, BB700, PerCP-Vio 700
	A	783/60		PE-Cy7	PE-A750, SBY800, PE-Vio 770
Red: 640	B	660/10		Alexa 647	Alexa 647, Far Red Live/Dead, SBR670, eF660, Vio Bright R667
	A	783/56		APC-Cy7	APC-H7, APC-eF780, SBR775, APC-Fire750, APC-Vio 770, L/D NIR, Zombie NIR